



Product Comparison Table



Photon
Bionova® Instant Auto-reader
 The quickest, newest and trendiest one



MiniBio
Bionova® Auto-reader
 SCBI incubation powerhouse at an affordable price



IC10/20FRLCD
Bionova® Auto-reader
 The most integrated and advanced auto-reader so far

NUMBER OF INCUBATION POSITIONS	2 for Self Contained Biological Indicators (SCBIs)	3 for Self Contained Biological Indicators (SCBIs)	12 for Self Contained Biological Indicators (SCBIs) 1 for Hygiene monitoring systems
COMPATIBLE INDICATORS	Bionova Photon® SCBIs	Bionova® SCBIs	Bionova® SCBIs Chemdye® Hygiene monitoring system indicators
COMPATIBLE MONITORED PROCESSES	STEAM Sterilization	STEAM (V)H2O2 EO FORM H2O2 Room Disinfection	STEAM (V)H2O2 EO FORM H2O2 Room Disinfection Hygiene monitoring system disinfection
AVAILABLE INCUBATIONS PROGRAMS	1 Available: Instant (7 seconds)	6 Available: Ultra Rapid (20 min) Super Rapid (30 min, 1 h) Rapid (2 hs, 3 hs, 4 hs)	7 Available: Ultra Rapid (20 min) Super Rapid (30 min, 1 h) Rapid (2 hs, 3 hs, 4 hs) PRO (15 min)
USER INTERFACE	Keypad Wireless assistant Mobile App	Keypad	LCD Touchscreen
TICKET FORMAT	Digital	Printed	Printed
CONNECTIVITY	USB Wi-Fi Bluetooth®	USB	Ethernet Wi-Fi (By using the Wi-link dongle)
ADVANCED SOFTWARE INTEGRATION	Fully compatible with Bionova® Classic and Cloud	Fully compatible with Bionova® Classic and Cloud	Fully compatible with Bionova® Classic & Cloud Auto-reader internal website accessible through internet browser