

Protein detection pen  
for cannulated devices

# Quantitative results in only 4 minutes with MiniPro!

Chemdye®  
PRO1 ENDO - SWE

The Chemdye® PRO1 ENDO Protein Detection Pen and Chemdye® SWE High-Absorption Swabs were specially designed to effectively verify the cleanliness of cannulated instruments, such as endoscopes, by detecting protein residues in their internal channels.



2.5 m swabs,  
available in 4  
diameters:



PRO1 ENDO is compatible with MINIPRO, IC10/20FR,  
IC10/20FRLCD auto-readers, and Bionova® Q software.

     
terragene.com

