

Sterilization efficiency
without concessions:

Safety guaranteed with Terragene PCDs




What is a PCD?

A **Process Challenge Device (PCD)** is designed to challenge the sterilization process under worst-case conditions, helping verify that the cycle performed effectively.



When the **PCD**
approves,
you know it's
safe.



Immediate results
at the end of the
sterilization cycle
— without
opening the
packages.



Why Terragene **PCDs**?



Regulated Precision
FDA Cleared



Error Reduction
Automated systems eliminate
human errors.



Versatility
Suitable for sterilizers of all sizes.



When should **PCDs** be used?

- Routine monitoring (daily / implant loads)
- Qualification testing after installation, repairs or suspected failures



Every load, every day, every test: Trust never takes a break.

Different PCD configurations.



Different sterilization processes require
different PCD configurations.

Contact our team to find the solution
that best fits your needs.



Maximize safety.
Optimize your process.
Trust Terragene.

Contact us today to learn
more and transform your
sterilization process.

terrogene.com/customer-services

