



CDWU-Z

Control of cleaning and washing of instruments and dental materials is highly important since its outcome affects the success of subsequent disinfection and/or sterilization processes.

- ✓ Reliable results in routine control of ultrasonic washers.
- ✓ High confidence and reliability to reveal deficiencies in the cavitation capacity of ultrasonic washing machines.
- ✓ Developed for releasing the ultrasonic washer together with the washing solution before the washing procedure.



UNPROCESSED



LOW CAVITATION
CAPACITY



HIGH CAVITATION
CAPACITY