

No textile bag



Digital traceability

ISO 11140-6:2022

Helix evolved

The next generation of Helix monitoring.

A Type 2 system with a reusable PCD for monitoring air removal and steam penetration in sterilization cycles involving narrow lumens.



Why upgrade to Helix Plus Kit?

- Protective reusable housing
- No textile bag
- Fewer handling steps
- Data Matrix identification
- Compatible with Trazanto Lens
- Integrated with Bionova® Q App

Compatible PCD-A indicators:

- PCD-A-3.5Y
- PCD-A-5.3Y
- PCD-A-7.0Y

Designed in conformity with:

- ISO 11140-1:2014
- ISO 11140-6:2022

